

GETTING STARTED

AC-3000 Series

Colposcope

Prepare for installation, learn the controls and organise the documentation for your AC-3000 Series.



Preparation guide · 14 September 2026

A preparation guide for trained healthcare teams. Use the operating instructions supplied for your unit and complete the appropriate training before clinical use.

01 Confirm the AC-3000 configuration

Check the stand, binocular, objective and magnification configuration on your order. The brochure lists standard and optional arrangements; not every option is supplied together.

02 Review the introduction

Watch the manufacturer video to identify the microscope and controls. Arrange supplier installation and a handover covering support, movement and the working position.

03 Learn the optical adjustments

Have the binocular, focus, magnification, illumination and tilting controls demonstrated on the delivered unit. Practise viewing a non-clinical target during training.

04 Set up documentation

Have the integrated camera, display and storage demonstrated. Request the operating instructions for your exact unit and camera version before clinical use.

Your Reference Material

[AC-3000 Series brochure](#)

Alltion Australia edition · October 2026 · PDF

[AC-3000 introduction](#)

Alltion

Brochures describe a range of configurations. Ask for the operating manual and setup instructions matched to your delivered unit.