

GETTING STARTED

AM-2000 Series

Dental Microscope

Prepare for installation, learn the controls and organise the documentation for your AM-2000 Series.



Preparation guide · 14 September 2026

A preparation guide for trained healthcare teams. Use the operating instructions supplied for your unit and complete the appropriate training before clinical use.

01 Match the delivered configuration

Check the model suffix, objective and camera against your order. The AM-2000 brochure covers several fixed-focus, variable-objective and camera configurations.

02 Review installation

Watch the series installation video and arrange the supplier handover for your mounting type. Floor, wall and ceiling installations need different setup arrangements.

03 Prepare your viewing position

Have the supplier demonstrate binocular adjustment, eyepieces, focus, magnification and illumination. Record the settings you need for a comfortable training session.

04 Check any imaging accessories

If a camera is fitted, identify its model and use its matching resources. Confirm display, storage and export with a test image before clinical use.

Your Reference Material

[AM-2000 Series brochure](#)

Alltion Australia edition · October 2026 · PDF

[AM-2000 installation introduction](#)

YouTube

[AM-2000 installation animation](#)

YouTube

[AM-2000 features and controls](#)

YouTube

Brochures describe a range of configurations. Ask for the operating manual and setup instructions matched to your delivered unit.

