

GETTING STARTED

AM-7400

Dental & Surgical Microscope

Prepare for the configuration review and supplier handover of your AM-7400.

Preparation guide · 20 October 2026



A preparation checklist for trained healthcare teams. It does not replace the operating instructions or training supplied with your unit.

- 01 Confirm the configuration**

Record whether the unit is Standard, Medium or High. Confirm zoom, illumination, eyepieces, objective and any additional accessories against your quotation.
- 02 Plan installation**

Confirm the mounting arrangement, available space, monitor and local power requirements with your supplier. Arrange assembly and balance checks by the appropriate personnel.
- 03 Review the fitted controls**

Ask the supplier to demonstrate focusing, zoom, filters, light-spot adjustment, touchscreen and arm-lock release on your delivered configuration.
- 04 Prepare documentation**

Request the model-specific operating instructions. During handover, review photo/video capture, storage and export using a non-clinical test.

Your Reference Material

[AM-7400 leaflet](#)

Alltion Australia edition · October 2026 · PDF

The leaflet covers Standard, Medium and High configurations and optional accessories. Confirm the exact supplied setup and current operating instructions with our team.