

GETTING STARTED

ANGEL 60

Dental & Surgical Microscope

Prepare for installation, configuration review and supplier training for your ANGEL 60.



Preparation guide · 20 October 2026

A preparation checklist for trained healthcare teams. It does not replace the operating instructions or training supplied with your unit.

- 01 Record your supplied setup**
Confirm the objective, eyepieces, imaging system, illumination options, stand and control accessories against the quotation and delivery record.
- 02 Arrange installation**
Ask the supplier to confirm the mounting position, balance-arm setup, clearance and power requirements. Obtain the correct assembly and operating instructions.
- 03 Learn the handgrip controls**
During handover, review focus, zoom, illumination, capture and electromagnetic-brake release. Confirm which optional functions and pedal controls are fitted.
- 04 Check imaging and support**
Use a non-clinical test to practise capture, storage and export. Keep the serial number, operating instructions and support contact with the equipment.

Your Reference Material

[ANGEL 60 leaflet](#)

Alltion Australia edition · October 2026 · PDF

The manufacturer leaflet describes standard and optional configurations. Confirm the delivered accessories, illumination functions and current operating instructions with our team.