



SOFTWARE & IMAGING

Colposcopy Software

+ HDMI Capture Card

Product code ALAHDMICAP

Professional colposcopy software with an HDMI capture card for the AC-2311 colposcope: image collection, atlas comparison, annotation, examination timing, reporting and case search.

AC-2311

Colposcope connection

175

Standard atlas images

11

Report types

Key features

- Set for the AC-2311: HDMI outside capture card, software and driver CD, HDMI cable and dongle.
- Displays, freezes, captures and deletes 1–99 colour images per inspection.
- Foot-switch image capture.
- Acetic-acid test timer.
- Atlas with 175 standard images; marking, remarks, measurement, zoom and area calculation.
- Clinical-symptom and treatment entry.
- Eleven report formats and RCI evaluation.
- Case creation, editing, deletion and search by number, name, age, date or fuzzy search.
- English, Spanish, Vietnamese and Chinese languages.

Specifications

PACKAGE & COMPATIBILITY

Product	HDMI capture card and professional colposcopy software set for the AC-2311 colposcope
Set contents	HDMI outside capture card; software and driver CD; HDMI cable; dongle — for connecting the AC-2311 to a desktop computer or laptop
Computer and monitor set	As above with 2 HDMI cables and an HDMI connector, for connecting the AC-2311 to a computer (or laptop) and a monitor at the same time

SOFTWARE & REPORTING

Image collection	Display, freeze, capture and delete; 1–99 colour images per inspection; foot-switch capture and roll browsing
Examination tools	Acetic-acid test timer; white balance and green filter
Atlas and annotation	175 standard atlas images; marking, remarks, measurement and partial or total magnification
Reports	11 kinds of reports; RCI evaluation
Case records	Add, delete, edit, review; search by case number, name, age or date; fuzzy search
Languages	English, Spanish, Vietnamese and Chinese
Software reference edition	v3.2.4 manual; features vary by software version

Specifications can change and may vary by configuration. Confirm the supplied specifications before ordering.

Alltion Australia · alltion.com.au · sales@alltion.com.au